



Jing Jing Pharmaceutical Co.,Ltd.

No. 88, Jingyi Road, Dacaozhuang Industrial Zone,

Ningjin County, Xingtai City, Hebei, P.R. China

Tel: +86-311-83078913 Fax: +86-311-83078913

E-mail: guolei@jjpharm.cn Website: <http://www.jjpharm.cn>



L- Alanine Specification Sheet

检测项目	标准
性状 Description	本品为白色结晶或结晶性粉末。 White crystalline or crystalline powder.
鉴别 Identification	样品的红外光吸收图谱与对照的图谱一致。 The infrared absorption spectrum of the sample is in agreement with that of the control .
PH	5.5 — 7.0
比旋度 (°) Optical rotation (°)	+13.7 — +15.1
氯化物 (%) Chloride (%)	≤0.05
硫酸盐 (%) Sulfate (%)	≤0.03
干燥失重 (%) Loss on drying (%)	≤0.2
炽灼残渣 (%) Residue on ignition (%)	≤0.15
铁 盐 (ppm) Iron (ppm)	≤30
含量测定 (%) Assay (%)	98.5-101.5